

# Experiment Brief

## Cipher system

### Title

Quantified backscattered electron analysis in the scanning electron microscope (SEM)

### Gatan Instrument Used

Cipher™ system, including the OnPoint™ backscattered electron detector and DigitalMicrograph® software.


### Background

In the scanning electron microscope (SEM), backscattered electron (BSE) imaging is used routinely to reveal the compositional variation within a specimen. The contrast in a BSE image results from the difference in average atomic number between analysis points. Typically, BSE imaging is performed qualitatively; however, quantitative BSE (qBSE) imaging may be performed under appropriate conditions to determine the mean atomic number ( $\bar{Z}$ ). Recently, interest in qBSE analysis has surged after researchers demonstrated the ability to reveal single-digit weight percentages of lithium in a sample quantitatively using qBSE and quantitative energy dispersive x-ray spectroscopy (EDS) analysis [1, 2]. In this article, we examine the use of the OnPoint BSE detector and DigitalMicrograph software (components of the Cipher system) for qBSE applications.


### Materials and Methods

qBSE analysis was performed on 55 high purity standards ( $4 < \bar{Z} < 83$ ). The samples were mechanically polished and coated with a 2.0 nm thick carbon layer using a PECS™ II system. Figure 1 shows the qBSE signal measured as a function of  $\bar{Z}$  (calculated using the modified electron approach after [3]). For compounds with  $\bar{Z} < 40$ , an excellent fit of the experimental data to the trendline was observed with few—if any—outliers.

A second sample was analyzed to investigate the accuracy of qBSE measurements in randomly orientated crystals. LiF particles, approximately cuboidal in shape and 1 – 15  $\mu\text{m}$  in size, were embedded in epoxy resin. The sample was milled to reveal the particle interior and coated with a 2.0 nm thick carbon layer using a PECS II system. qBSE analysis of the particles was performed, Figure 2. The  $\bar{Z}$  determined experimentally was  $7.15 \pm 0.19$ , comparing favorably with the calculated value of 7.10. Analysis of the lithium content by Cipher reveals the mean lithium content to be  $48.8 \pm 0.04$  at. % and demonstrates that channeling contrast is not a significant factor in determining the lithium content under these conditions.



**Figure 1.** The plot of normalized backscattered electron grey levels against mean atomic number was calculated using the modified electron approach (after [3]). The circles are experimental data collected at an accelerating voltage of 15 kV, and the dotted line is the exponential fit function.



**Figure 2.** Mean atomic number of 250 LiF particles as determined by qBSE analysis using the Cipher system. Collected at an accelerating voltage of 15 kV.

### Summary

Precise determination of a sample's mean atomic number ( $\bar{Z}$ ) by quantitative backscattered electron analysis was demonstrated for a wide range of materials of  $\bar{Z} < 40$ . Crystal orientation was found to have minimal effect on qBSE results enabling the lithium content of LiF particles to be determined with an accuracy of 1.2 at. % (<1 wt. %) by the Cipher system. These results promise exciting characterization possibilities of lithiated materials.

**Gatan, Inc.** is the world's leading manufacturer of instrumentation and software used to enhance and extend electron microscopes—from specimen preparation and manipulation to imaging and analysis.

### Credit(s)

- [1] J.A. Österreicher et al., Scripta Materialia 194 (2021), 113664
- [2] J. Lee et al., Microsc. Microanal. 28 (2022), p548-550
- [3] J. Donovan et al., Microsc. Microanal. 9 (2003), p202