

# Experiment Brief

## Metro counting camera

### Title

Enhancing MicroED/3DED analysis with direct detection electron counting cameras

### Gatan instrument used


The Metro® camera enables low-dose imaging and **diffraction** via **real-time electron counting** with a simple user interface.

### Background


PPEA, a semiconducting molecule, exhibits a high affinity to graphene surfaces. The PPEA molecule forms very small, needle-like crystals that measure several micrometers in length and have a width ranging from 100 to 200 nm. Due to their size, these crystals are too small for analysis using single-crystal x-ray structure techniques. Moreover, the beam-sensitive nature of PPEA requires low-dose techniques in the transmission electron microscope (TEM) to avoid damaging the crystallites, making continuous diffraction tomography (MicroED/3DED) the method of choice for structure determination.

### Materials and methods

PPEA crystals were grown and prepared for TEM using the methods described in Reference 1. Data was acquired at 200 kV with a Metro camera. A continuous tilt series from  $-30^\circ$  to  $30^\circ$  was collected using a tilt speed of  $\sim 1^\circ/\text{s}$  tilt and a dose rate of  $0.0311 \text{ e}^-/\text{\AA}^2/\text{s}$  to ensure a total dose of  $< 2 \text{ e}^-/\text{\AA}^2$  for the tilt series acquisition.



**Figure A:** Molecular structure and **B:** Conformation within the crystal structure of the 9,10-bis-((perchloro-phenyl)-ethynyl)anthracene (PPEA). **Figure C:** Single frame from a continuous diffraction tomography dataset showcasing higher order reflections with resolutions exceeding 1 Å. The full dataset video is available at [https://youtu.be/Nr3v4vTa\\_IY](https://youtu.be/Nr3v4vTa_IY).



### Summary

This study highlights the effectiveness of the Metro electron counting camera for MicroED/3DED experiments without needing a beam stop. A key achievement here is maintaining a low electron dose rate during full tomogram data collection, which is crucial for studying sensitive materials such as PPEA without altering their natural state. Further work could leverage the Latitude® D software for high-throughput automated data collection of continuous diffraction tomography datasets.

### Reference

1) Gorelik, T.E., Ulmer, A., Schleper, A.L., Kuehne, A.J.C., Crystal Structure of 9,10-bis-((perchloro-phenyl)-ethynyl)anthracene Determined from Three-Dimensional Electron Diffraction Data, Z. Kristallogr. (2023), <https://doi.org/10.1515/zkri-2023-0009>

### Credit(s)

Special thanks to Dr. Tatiana Gorelik at the Helmholtz Centre for Infection Research.

